

Blepharitis

Blephaclean Wipes

Blephasol



FREE FROM PRESERVATIVES & PERFUME

FREE FROM PRESERVATIVES & PERFUME




Blephasol and Blephaclean Wipes, for effective control and treatment of blepharitis and eye lid hygiene.

Blephasol and Blephaclean Wipes, how they work.

Lipid secretions adsorbed into the centre of the micelle

Hydrophilic pole

Hydrophobic centre



Thanks to the bipolar structure of micelles, micelle solutions eliminate lipids by a simple mechanical phenomenon: ADSORPTION

OSD Management

To enhance the effectiveness of tear lubricants⁽¹⁾

Surgical Prophylaxis

To prevent pre and post-surgery ocular infection⁽²⁾

(1) - Benitez del Castillo JM, Arangué C, García Sánchez J. Corneal epithelial permeability and dry eye treatment. Third international conference on the lacrimal gland tear film and dry eye syndromes : basic science and clinical relevance. Conference program and abstract book. Cornea 2001, 19, suppl 6, S76.

(2) - Hueso Abancens JR, Mengual Verdu E, Schargel Palacios K, García Sánchez J, Placeres Daban J, García Conca V, Zorraquino Martí A. Modification of the conjunctival flora caused by eyelid cleaning solutions.

Formulated by Dermatologists • Tested by Ophthalmologists

SpectrumThea
Fernbank House Springwood Way
Macclesfield Cheshire SK10 2XA
Tel: 0845 521 1290 Fax: 0845 521 1290

www.spectrumthea.co.uk

